



## How To Use This Catalog Page Reference

**Product Type** ————— *Fixed Curve Diagnostic Catheters*

**French Size and Number of Electrodes** ————— **5F WEBSTER® Quadrapolar Catheter with Redel Connector**

**Electrode Spacing** ————— 

**Type of Connector & Associated Cable** ————— **Redel 10-Pin Connector** Associated Cable:  
Catalog No: C6MRMST4SA  
European Order No: 39F-55R

**New U.S. Catalog Number** ————— 


New Catalog No.	Made to Order	Old Catalog No.	Indication	Tip Electrode	Size	Curve/Color	Spacing	Insertion Length	European Order No.
F5QA002RT	•	F5-QA-002-RT	RA, RV	1 mm	5F	A/Yellow	2 mm	110 cm	37C-00R
F5QA252RT		F5-QA-252-RT	RA, RV, His	1 mm	5F	A/Yellow	2-5-2 mm	110 cm	37C-08R
F5QA005RT		F5-QA-005-RT	RA, RV	1 mm	5F	A/Yellow	5 mm	110 cm	37C-03R
F5QA010RT	•	F5-QA-010-RT	RA, RV, His	1 mm	5F	A/Yellow	10 mm	110 cm	37C-05R
F5QD002RT	•	F5-QD-002-RT	RV	1 mm	5F	D/Blue	2 mm	110 cm	37C-30R
F5QD252RT	•	F5-QD-252-RT	RA, RV	1 mm	5F	D/Blue	2-5-2 mm	110 cm	37C-38R
F5QD005RT		F5-QD-005-RT	RA, RV	1 mm	5F	D/Blue	5 mm	110 cm	37C-33R
F5QD010RT	•	F5-QD-010-RT	RV, His	1 mm	5F	D/Blue	10 mm	110 cm	37C-35R
F5QF002RT	•	F5-QF-002-RT	RA, RV	1 mm	5F	F/Black	2 mm	110 cm	37C-50R
F5QF252RT	•	F5-QF-252-RT	RA, RV, His	1 mm	5F	F/Black	2-5-2 mm	110 cm	37C-58R
F5QF005RT		F5-QF-005-RT	RA, RV	1 mm	5F	F/Black	5 mm	110 cm	37C-53R
F5QF010RT		F5-QF-010-RT	RA, RV, His	1 mm	5F	F/Black	10 mm	110 cm	37C-55R
F5QL005RT	•	F5-QL-005-RT	RV, His	1 mm	5F	L	5 mm	110 cm	37C-83R

**Made to Order Status** —————

**Old U.S. Catalog Number** —————

**European Order Number** —————

Customers may expect orders to ship the same day as placed. Product availability will be confirmed at time of order placement.

 **Biosense Webster**  
a Johnson & Johnson company

Your Biosense Webster Inc. Sales Representative is available to assist you by answering questions about our products and service.