

CELSIUS® DS Temperature Sensing Ablation Catheters 7F Quadrapolar with 8mm Tip/Thermocouple



Redel 10-Pin Connector

Associated Cables:
Catalog No: C10MRMSTKDTCS (Stockert 70 Generator)
European Order No: 39E-68R



New Catalog No.	Made to Order	Old Catalog No.	Tip Electrode	Size	Curve/Color	Spacing	Insertion Length	European Order No.
D7TCB8L162RT		D7TC-B8L-162-RT	8 mm	7F	B-Type/Red	1-6-2 mm	115 cm	35S-17R
D7TCD8L162RT		D7TC-D8L-162-RT	8 mm	7F	D-Type/Blue	1-6-2 mm	115 cm	35S-37R
D7TCF8L162RT		D7TC-F8L-162-RT	8 mm	7F	F-Type/Orange	1-6-2 mm	115 cm	35S-57R
D7TCJ8L162RT		D7TC-J8L-162-RT	8 mm	7F	J-Type/Black	1-6-2 mm	115 cm	N/A