

CELSIUS® Temperature Sensing Ablation Catheters 6F and 7F Quadrapolar with Non-Braided Tip/Thermistor



Redel 10-Pin Connector




Associated Cables:
Catalog No: C10MR10MSTKS (Stockert 70 Generator)
European Order No: 39E-43R



Catalog No: C6MR10EPTST (EPT-1000* TC Generator)
European Order No: 39E-32R



Catalog No: C6MR10EPTRS (EPT-1000* TC Generator)
European Order No: 39E-69R



New Catalog No.	Made to Order	Old Catalog No.	Tip Electrode	Size	Curve/Color	Spacing	Insertion Length	European Order No.
D6TAG252RT	●	D6T-AG-252-RT	4 mm Grooved	6F	A-Type/Yellow	2-5-2 mm	92 cm	N/A
D6TAL252RT		D6T-AL-252-RT	4 mm	6F	A-Type/Yellow	2-5-2 mm	92 cm	35H-07R
D6TBL252RT	●	D6T-BL-252-RT	4 mm	6F	B-Type/Red	2-5-2 mm	92 cm	35H-17R
D6TDG252RT	●	D6T-DG-252-RT	4 mm Grooved	6F	D-Type/Blue	2-5-2 mm	92 cm	N/A
D6TDL252RT		D6T-DL-252-RT	4 mm	6F	D-Type/Blue	2-5-2 mm	92 cm	35H-37R
D7TAL252RT		D7T-AL-252-RT	4 mm	7F	A-Type/Yellow	2-5-2 mm	115 cm	360-07R
D7TBL252RT		D7T-BL-252-RT	4 mm	7F	B-Type/Red	2-5-2 mm	115 cm	360-17R
D7TCG252RT	●	D7T-CG-252-RT	4 mm Grooved	7F	C-Type/Green	2-5-2 mm	115 cm	N/A
D7TCL252RT		D7T-CL-252-RT	4 mm	7F	C-Type/Green	2-5-2 mm	115 cm	360-27R
D7TDG252RT		D7T-DG-252-RT	4 mm Grooved	7F	D-Type/Blue	2-5-2 mm	115 cm	N/A
D7TDL252RT		D7T-DL-252-RT	4 mm	7F	D-Type/Blue	2-5-2 mm	115 cm	360-37R
D7TEG252RT	●	D7T-EG-252-RT	4 mm Grooved	7F	E-Type/White	2-5-2 mm	115 cm	N/A
D7TEL252RT		D7T-EL-252-RT	4 mm	7F	E-Type/White	2-5-2 mm	115 cm	360-47R
D7TFG252RT	●	D7T-FG-252-RT	4 mm Grooved	7F	F-Type/Orange	2-5-2 mm	115 cm	N/A
D7TFL252RT		D7T-FL-252-RT	4 mm	7F	F-Type/Orange	2-5-2 mm	115 cm	360-57R