

NAVI^{STAR}® Soft Tip Temperature Sensing Navigation and Ablation Catheter 7F Quadrapolar with 4 mm Tip/Thermocouple



Hypertronics 25-Pin Connector

Associated Cables: CARTO™ System Junction Box
Catalog No: C5MHNAMHS (Red)
European Order No: 39E-58M



CARTO™ System Junction Box to Generator
Catalog No: C10MR10MSTKS (Stockert 70 Generator)
European Order No: 39E-43R



Catalog No: C6MR10MLTCS (Medtronic ATAKR® Generator)
European Order No: 39E-23R



New Catalog No.	Made to Order	Old Catalog No.	Tip Electrode	Size	Curve/Color	Spacing	Insertion Length	European Order No.
118334S	●	1183-34-S	4 mm	7F	B-Type/Red	1-7-4 mm	115 cm	SW1183-035
118330S	●	1183-30-S	4 mm	7F	D-Type/Blue	1-7-4 mm	115 cm	SW1183-030
118336S	●	1183-36-S	4 mm	7F	E-Type/White	1-7-4 mm	115 cm	SW1183-036
118331S	●	1183-31-S	4 mm	7F	F-Type/Orange	1-7-4 mm	115 cm	SW1183-031
118332S	●	1183-32-S	4 mm	7F	J-Type/Black	1-7-4 mm	115 cm	SW1183-032