

## NAVISTAR® DS Temperature Sensing Navigation and Ablation Catheters 7F Quadrapolar with 8 mm Tip/Thermocouple



### Hypertronics 25-Pin Connector

Associated Cables: CARTO™ System Junction Box  
Catalog No: C5MHDTMHS (Green)  
European Order No: 39F-12M



CARTO™ System Junction Box to Generator  
Catalog No: C10MRMSTKDTCS (Stockert 70 Generator)  
European Order No: 39E-68R



New Catalog No.	Made to Order	Old Catalog No.	Tip Electrode	Size	Curve/Color	Spacing	Insertion Length	European Order No.
NS7TCB8L174HS		NS7TC-B8L-174-HS	8 mm	7F	B Type/Red	1-7-4 mm	115 cm	34J-17M
NS7TCC8L174HS	●	NS7TC-C8L-174-HS	8 mm	8/7F	C Type/Green	1-7-4 mm	115 cm	34J-27M
NS7TCD8L174HS		NS7TC-D8L-174-HS	8 mm	7F	D Type/Blue	1-7-4 mm	115 cm	34J-37M
NS7TCF8L174HS		NS7TC-F8L-174-HS	8 mm	7F	F Type/Orange	1-7-4 mm	115 cm	34J-57M
NS7TCJ8L174HS		NS7TC-J8L-174-HS	8 mm	7F	J Type/Black	1-7-4 mm	115 cm	34J-J7M