

QwikStar® DS Temperature Sensing Navigation and Ablation Catheters 7F 26-Pole with 8 mm Tip/Thermocouple



Microhelix 70-Pin Connector

Associated Cable: CARTO™ XP System PIU

Catalog No: C5MYQFM

European Order No: C5MYQFM



CARTO™ XP System PIU to Generator

Catalog No: C10MR10MSTKS (Stockert 70 Generator)

European Order No: 39E-43R



Catalog No: C6MR10MLTCS (Medtronic ATAKR® Generator)

European Order No: 39E-23R



Catalog No.	Made to Order	Tip Electrode	Size	Curve/Color	Spacing	Insertion Length	European Order No.
QS7TCB8M		8 mm	7F	B-Type/Red	1-9-9-9-9-9 mm	115 cm	QS7TCB8M
QS7TCD8M		8 mm	7F	D-Type/Blue	1-9-9-9-9-9 mm	115 cm	QS7TCD8M
QS7TCF8M		8 mm	7F	F-Type/Orange	1-9-9-9-9-9 mm	115 cm	QS7TCF8M
QS7TCJ8M		8 mm	7F	J-Type/Black	1-9-9-9-9-9 mm	115 cm	QS7TCJ8M