

BIOSENSE WEBSTER CABLE QUICK REFERENCE

Fixed and Deflectable Diagnostic Catheters		
Catheter Connector Redel 10-Pin	Catheter Type Fixed Catheters and Deflectable Catheters with product codes ending in RT	
Catalog No. C6MRMST4SA	Catheter Mating Connector Redel 10-Pin	Recorder/Wiring Connector Four (4) Shielded Tip Pins
Catalog No. C6MRMST6SA	Catheter Mating Connector Redel 10-Pin	Recorder/Wiring Connector Six (6) Shielded Tip Pins
Catalog No. C6MRMST8SA	Catheter Mating Connector Redel 10-Pin	Recorder/Wiring Connector Eight (8) Shielded Tip Pins
Catalog No. C6MRMST10SA	Catheter Mating Connector Redel 10-Pin	Recorder/Wiring Connector Ten (10) Shielded Tip Pins

Non-Navigation Therapeutic Catheters		
Catheter Connector Redel 10-Pin	Catheter Type Catheters with product codes ending in RT	
Catalog No. C6MR10MSTKS C10MR10MSTKS Generator Model Stockert 70 (Thermocouple or Thermistor)	Catheter Mating Connector Redel 10-Pin	RF Generator Connector Redel 10-Pin
Catalog No. C6MR10EPTST Generator Model EPT-1000 (Thermistor)	Catheter Mating Connector Redel 10-Pin	RF Generator Connector Shielded Tip Pins
Catalog No. C6MR10EPTRS Generator Model EPT-1000 (Thermistor)	Catheter Mating Connector Redel 10-Pin	RF Generator Connector Redel 9-Pin
Catalog No. C6MR10MLTCS Generator Model Medtronic Atakr (Thermocouple)	Catheter Mating Connector Redel 10-Pin	RF Generator Connector Lemo 10-Pin
Catheter Connector Redel 10-Pin	Catheter Type DS Catheters with product codes ending in RT	
Catalog No. C6MRMSTKDTCS C10MRMSTKDTCS Generator Model Stockert 70 (Dual Thermocouples)	Catheter Mating Connector Redel 10-Pin	RF Generator Connector Redel 10-Pin

Navigation Therapeutic Catheters		
Catheter Connector Hypertronics 25-Pin	Catheter Type NAVISTAR® DS Catheters with product codes ending in HS	
Catalog No. C5MHDTCMHS	Catheter Mating Connector Hypertronics 25-Pin	CARTO™ System Junction Box Connector Hypertronics 25-Pin
Catheter Connector Hypertronics 25-Pin	Catheter Type NAVISTAR® and REFSTAR™ Catheters with product codes ending in HS	
Catalog No. C5MHNAMVMS*	Catheter Mating Connector Hypertronics 25-Pin	CARTO™ System Junction Box Connector Hypertronics 25-Pin
Catalog No. C5MHREFMHS*	Catheter Mating Connector Hypertronics 25-Pin	CARTO™ System Junction Box Connector Hypertronics 25-Pin
Catheter Connector Microhelix 70-Pin	Catheter Type QWIKSTAR® Catheters with product codes ending in M	
Catalog No. C5MYQFM	Catheter Mating Connector Microhelix 70-Pin	CARTO™ System Junction Box Connector Hypertronics 70-Pin

CARTO™ System Junction Box		
Catalog No. C6MR10MSTKS C10MR10MSTKS Generator Model Stockert 70 (Thermocouple or Thermistor)	CARTO™ System Junction Box Connector Redel 10-Pin	RF Generator Connector Redel 10-Pin
Catalog No. C6MRMSTKDTCS C10MRMSTKDTCS Generator Model Stockert 70 (Dual Thermocouples)	CARTO™ System Junction Box Connector Redel 10-Pin	RF Generator Connector Redel 10-Pin
Catalog No. C6MR10EPTST Generator Model EPT-1000 (Thermistor)	CARTO™ System Junction Box Connector Redel 10-Pin	RF Generator Connector Shielded Tip Pins
Catalog No. C6MR10EPTRS Generator Model EPT-1000 (Thermistor)	CARTO™ System Junction Box Connector Redel 10-Pin	RF Generator Connector Redel 9-Pin
Catalog No. C6MR10MLTCS Generator Model Medtronic Atakr (Thermocouple)	CARTO™ System Junction Box Connector Redel 10-Pin	RF Generator Connector Lemo 10-Pin

*Cables are electrically identical and can be interchanged.

CABLES

Catalog No.	Length	Catheter Connector Type	Electrogram Conductors	Equipment Connector Type	Equipment Interface
C6MR10MSTKS	6 Feet	Redel 10-Pin	4	Redel 10-Pin	Stockert 70 Generator
C10MR10MSTKS	10 Feet	Redel 10-Pin	4	Redel 10-Pin	Stockert 70 Generator
C6MRMSTKDTCS	6 Feet	Redel 10-Pin	4	Redel 10-Pin	Stockert 70 Generator
C10MRMSTKDTCS	10 Feet	Redel 10-Pin	4	Redel 10-Pin	Stockert 70 Generator
C6MR10EPTST	6 Feet	Redel 10-Pin	4	Shielded Tip Pins	EPT-1000 Generator
C6MR10EPTRS	6 Feet	Redel 10-Pin	4	Redel 9-Pin	EPT-1000 Generator
C6MR10MLTCS	6 Feet	Redel 10-Pin	4	Lemo 10-Pin	Medtronic Atakr Generator
C5MHNAVMHS	5 Feet	Hypertronics 25-Pin	4	Hypertronics 25-Pin	CARTO™ System Junction Box
C5MHREFMHS	5 Feet	Hypertronics 25-Pin	4	Hypertronics 25-Pin	CARTO™ System Junction Box
C5MHDTCMHS	5 Feet	Hypertronics 25-Pin	4	Hypertronics 25-Pin	CARTO™ System Junction Box
C5MYQFM	5 Feet	Microhelix 70-Pin	4	Hypertronics 70-Pin	CARTO™ System Junction Box
C6MRMST4SA	● 6 Feet	Redel 10-Pin	4	Shielded Tip Pins	
C6MRMST6SA	● 6 Feet	Redel 10-Pin	6	Shielded Tip Pins	
C6MRMST8SA	● 6 Feet	Redel 10-Pin	8	Shielded Tip Pins	
C6MRMST10SA	● 6 Feet	Redel 10-Pin	10	Shielded Tip Pins	

● Autoclavable

Caution: Federal law restricts these devices to sale by or on the order of a physician.
Please refer to the complete product information accompanying each device.

The third-party trademarks used herein are trademarks of their respective owners.

©2006 Biosense Webster, Inc. All rights reserved. Order No. **S0049** Printed in USA 000000000.0



Biosense Webster, Inc.

3333 Diamond Canyon Rd.

Diamond Bar, CA 91765, USA

Tel: 909-839-8500

Tel: 800-729-9010

Fax: 909-468-2905

www.biosensewebster.com