

Cables for CARTO™ System Junction Box/Patient Interface Unit


Connects to Junction Box



Connects to Generator



New Catalog No.	Made to Order	Old Catalog No.	Length	Junction Box Connector Type	Electrogram Conductors	Generator Connector Type	Typical Usage	European Order No.
C10MRMSTKDTCS		C10-MR10/MSTKDTC-S	10 Feet	Redel 10-pin	4	Redel 10-pin	Stockert 70 Generator	39E-68R
C10MR10MSTKS		C10-MR10/MSTK-S	10 Feet	Redel 10-pin	4	Redel 10-pin	Stockert 70 Generator	39E-43R
C6MR10EPTRS		C6-MR10/EPTR-S	6 Feet	Redel 10-pin	4	Redel 9-pin	EPT-1000* TC Generator	39E-69R



New Catalog No.	Made to Order	Old Catalog No.	Length	Junction Box Connector Type	Electrogram Conductors	Generator Connector Type	Typical Usage	European Order No.
C6MR10EPTST	●	C6-MR10/EPT-ST	6 Feet	Redel 10-pin	4	Shielded Tip Pins	EPT-1000* TC Generator	39E-32R



New Catalog No.	Made to Order	Old Catalog No.	Length	Junction Box Connector Type	Electrogram Conductors	Generator Connector Type	Typical Usage	European Order No.
C6MR10MLTCS		C6-MR10/MLTC-S	6 Feet	Redel 10-pin	4	Lemo 10-pin	Medtronic ATAKR* Generator	39E-23R

Connects To Catheter



Connects To Junction Box



New Catalog No.	Made to Order	Old Catalog No.	Length	Junction Box Connector Type	Electrogram Conductors	Generator Connector Type	Typical Usage	European Order No.
C6MR10EPTRS		C6-MR10/EPTR-S	6 Feet	Redel 10-pin	4	Redel 9-pin	EPT-1000* TC Generator	39E-69R