


## Cables for Therapeutic Catheters with Redel Connector (RT)\*




Connects To Catheter



Connects to Generator



New Catalog No.	Made to Order	Old Catalog No.	Length	Catheter Connector Type	Electrogram Conductors	Generator Connector Type	Typical Usage	European Order No.
<b>C10MR10MSTKS</b>		<b>C10-MR10/MSTK-S</b>	10 Feet	Redel 10-Pin	4	Redel 10-Pin	Stockert 70 Generator	39E-43R
<b>C6MR10EPTRS</b>		<b>C6-MR10/EPTR-S</b>	6 Feet	Redel 10-Pin	4	Redel 9-Pin	EPT-1000** TC Generator	39E-69R



New Catalog No.	Made to Order	Old Catalog No.	Length	Catheter Connector Type	Electrogram Conductors	Generator Connector Type	Typical Usage	European Order No.
<b>C6MR10EPTST</b>	●	<b>C6-MR10/EPT-ST</b>	6 Feet	Redel 10-Pin	4	Shielded Tip Pins	EPT-1000** TC Generator	39E-32R



New Catalog No.	Made to Order	Old Catalog No.	Length	Catheter Connector Type	Electrogram Conductors	Generator Connector Type	Typical Usage	European Order No.
<b>C6MR10MLTCS</b>		<b>C6-MR10/MLTC-S</b>	6 Feet	Redel 10-Pin	4	Lemo 10-Pin	Medtronic ATAKR** Generator	39E-23R