


Cables for NAVI STAR[®] Catheters with Hypertronics Connector (HS)*



Connects To Catheter

Connects To Junction Box



New Catalog No.	Made to Order	Old Catalog No.	Length	Catheter Connector Type	Electrogram Conductors	Junction Box Connector Type	Typical Usage	European Order No.
C5MHNAVMHS		C5-MH/NAVMH-S	5 Feet	Hypertronics 25-Pin (Red)	4	Hypertronics 25-Pin (Red)	w/CARTO™ EP Navigation System Junction Box	39E-58M